




As a Full Service Digital Dental Laboratory, Centreline Dental Lab accepts files from all scanners. Find your scanner below, and follow the instructions, in order to start sending us your cases today!

MEDIT Link

- Go to www.meditLink.com
- Log in with your username and password
- Click on Partners search info@cldentallab.com
- Send a request by clicking 

Once we accept your request you will be able to “order” your case after scanning

iTero

- Call iTero support at 1.800.577.8767
- Select option 1, then option 1
- Request to have our lab added to scanner
- Provide our lab name, [Centreline Dental Lab](mailto:info@cldentallab.com), and the iTero customer number 889611
- On the scanner - go to settings and select “Sync Configuration”

Carestream

- Visit Carestream connect on your scanner
- Search for [Centreline Dental Lab \(info@cldentallab.com\)](mailto:info@cldentallab.com)
- Add Lab
- Select [Centreline Dental Lab](mailto:info@cldentallab.com) when submitting scans or raw .STL files
- Email your favorite digital lab, [Centreline Dental Lab](mailto:info@cldentallab.com)



- Log in to Cerec Connect at www.sirona-connect.com
- Select “My Cerec Connect” then “Edit Account”
- Under “My Account” select “My Favorite Laboratories”
- Type in name [Centreline Dental Lab \(info@cldentallab.com\)](mailto:info@cldentallab.com) in the field for “Company name” then the Postal Code
- Click “Find”
- Check the box beside [Centreline Dental Lab](mailto:info@cldentallab.com)
- Click “Add” to complete the registration



- Useful link: [3Shape TRIOS - Making connections - YouTube](#)
 - Log in to 3Shape software or link: <https://us.3shapecommunicate.com>
 - Click “More” and select “Settings”
 - Then in the menu click “Connections”
 - Click on “Labs” then the “Add”
 - Search [Shanto Dental Ceramics](mailto:info@shantodental.com) (our sister lab that supports us) by email address: info@shantodental.com
 - Click “Connect” to send a connection request
 - We will approve your request ASAP
- Once we approved your connection request it will be added to your list of connected Labs
- Click on “Patients” and “Save” your changes

